



# PORTABLE ORAL IRRIGATOR

## Product Manual

CF500

### 4 Modes Adjustment



#### Soft Mode

Sensitive or first time user



#### Strong Mode

Meet daily cleaning



#### Pulse Mode

Sensitive or first time user




#### Jog Mode

Easy to control strength



## Product structure diagram



- |                     |                                            |                         |
|---------------------|--------------------------------------------|-------------------------|
| ① Jet Nozzle        | ② Charging port                            | ③ Machine main body     |
| ④ Water-tank        | ⑤ Water inlet                              | ⑥ Nozzle release button |
| ⑦ Switch button     | ⑧ Mode switching                           | ⑨ Indicator lamp        |
| ⑩ Soft Mode         | ⑪ Strong Mode                              | ⑫ Pulse mode            |
| ⑬ Point-motion mode | ⑭ Charging indication / working indication |                         |

## Warm Tips

It is recommended to use it in pulse mode or gentle mode. After getting used to it, spray it in strong mode. When using it, do not leave a long distance and wash away food debris close to the gap between teeth. At the beginning, use with the mouth closed, not only won't splash water everywhere, but also because there is water in the mouth, the gums won't hurt.



## Safety instructions

### Dangerous:

- ① If you do not understand the operating procedures and safety instructions, please do not worry to use the product to avoid improper operation to damage.
- ② Do not open the product or change the parts of the product, so as not to damage the product and void the warranty.
- ③ Do not use water above 50 °C into the water tank and clean the mouth to avoid scald.
- ④ If the power cord or plug is damaged, replace it before use.
- ⑤ Do not put the product into fire or heat, nor charge, use, or place it at high temperatures.

### Warning:

- ① Do not use this product for any purpose other than oral cleansing.
- ② Do not press the nozzle against the gums or teeth too hard, so as not to damage the gums.
- ③ Do not immerse the product in water for a long time, so as not to fail or cause fire, electric shock and other dangers.
- ④ Do not damage, destroy, modify the power cord, do not stretch, bend, twist the power cord, do not press heavy objects on the power cord or put the power cord between objects, and do not let the product and the power cord close to hot objects.
- ⑤ Do not store in hot or humid places, and avoid direct sunlight.
- ⑥ Do not charge the product in a hot or humid place, such as the bathroom, toilet.
- ⑦ If not used for a long time, please unplug the power plug and charging plug, so as not to cause fire or danger.

### Notice:

- ① This product has no replaceable or repairable parts. If you need to repair, please return it to the manufacturer.
- ② Do not use it if there is no nozzle.
- ③ This product is not suitable for children, infants, etc. without operators, unless they receive supervision or guidance by their safety supervisor in the use of this product. Please place outside a child or baby, not used as a toy.
- ④ Patients with periodontal disease or have done oral surgery in nearly 2 months, please consult your doctor before use.



## Technical parameters


|                  |                         |                 |                 |
|------------------|-------------------------|-----------------|-----------------|
| Product Name     | Portable oral irrigator | Levels          | 4 Levels        |
| Nozzle diameter  | 0.8mm                   | Working current | 1.35A           |
| Waterproof grade | IPX7                    | Battery type    | lithium battery |
| Input            | DC5.0V                  | Rated power     | 5.5W            |
| Charging port    | TYPE-C                  | Charging time   | About 3 hours   |
| Product Size     | 68*45*148mm             | Tank capacity   | 160mL           |

## Fault diagnosis


| Problem                                          | Possible causes                                                      | Solutions                                      |
|--------------------------------------------------|----------------------------------------------------------------------|------------------------------------------------|
| Not working                                      | New purchase or more months not used                                 | Please charge it for more than 2 hours         |
| Low water pressure at working                    | There is something inside the nozzle blocking the nozzle deformation | Replace the nozzle                             |
| It only lasts a short time when fully charged    | Battery comes to the end of its service life                         | Please charge for more than 2 hours before use |
|                                                  | Insufficient charging time                                           |                                                |
| Charging indicator light is not on when charging | Power plug is not plugged in                                         | Please connect the power supply properly       |
|                                                  | Product charging circuit is faulty                                   |                                                |

# Product usage method


1. Install the nozzle.




2. Open the water tank, open the water injection port and add water.



3. Put the nozzle into the mouth and press the switch, start to use.




4. Regular cleaning: pull out the water tank box and clean the water tank.





## Fear Description

| Working mode                          | Mode switch                                                                                   | Mode indicator light colour | Applicable people                                 |
|---------------------------------------|-----------------------------------------------------------------------------------------------|-----------------------------|---------------------------------------------------|
| Soft Mode                             | Click the mode key to switch different gears, and then press the open key to work             | White、Green                 | Suitable for sensitive teeth and initial users    |
| Strong Mode                           | Click the mode key to switch different gears, and then press the open key to work             | White、Green                 | Meet most People's Daily cleaning                 |
| Pulse mode                            | Click the mode key to switch different gears, and then press the open key to work             | White、Green                 | Tooth sensitivity, oral sensitivity               |
| Point-motion mode                     | Click the mode key to switch different gears, and then press the open key to work             | White、Green                 | Convenient and comfortable, good strength control |
| Charge indicator<br>Work instructions | Respiratory flashing when charging, full of long light, low voltage reminder (slow flash 10S) |                             |                                                   |



### ① Jet Nozzle:

Press the nozzle unplug button, remove the nozzle, clean with water, and dry with a rag.

- Press the nozzle unplug button, remove the nozzle, clean with water, and dry with a rag.
- If the tip of the nozzle is deformed, replace a new nozzle or a fault can result.

### ② Machine main body:

Remove the dirt with a cloth dipped in liquid soap, and then wipe it clean with a dry cloth. If the body is dirty, wipe it with water.

- The dentist recommends replacing the nozzles every six months.
- If the tip of the nozzle is deformed, replace a new nozzle or a fault can result.

### ③ Water filling nozzle:

Wipe off the dirt or it will cause a leak.

### ④ Water-tank:

Unplug the tank box regularly and clean the tank.

- Wash off the internal dirt with items like brushes and wipe it off with a dry rag.
- Water tank shall be dried for long service.



## Environmental protection and resource recovery

When the product reaches its service life, it may be necessary to remove the internal battery for recycling and proper disposal according to the laws and regulations of the country or region where it is located.

If you cannot remove the battery, please submit the product to the government-designated recycling center.

This product does not support users to replace the battery by themselves. It will damage the function of the product during the forced battery removal. Please operate carefully.

### Built-in rechargeable battery removal method

This product contains rechargeable batteries. Please remove the batteries from the product before consumers decide to discard the product.

1. Decompose the product and turn off the power before taking out the battery.
2. Remove the battery step by step according to the sequence steps ① — ④ as shown in the figure on the right.
3. Batteries must be discarded in places or recycling units specified by the government.

